



Brewer
more possibilities

ASSIST POWER PROCEDURE TABLE

BrewerAssist

POWER PROCEDURE TABLE Infinitely adjustable for more possibilities.

The Brewer Assist Power Procedure Table is clearly the:

- most technologically advanced
- lowest for lateral patient transfer
- greatest weight bearing
- most adaptable
- most extendable
- most user-friendly
- most patient-friendly
- best warranted table in the industry.

**See why Brewer means...
more possibilities.**





Easy wheelchair patient transfer at 19"



Supports 450 lbs.



Exclusive, easy-to-adjust three-point pivot headrest system



Convenient table paper storage



Comfortable, ergonomic upholstery; easy to clean surface



Footrest extends an additional 5"



Standard



Optional

Intuitive foot and handheld controls

INFINITELY ADJUSTABLE



Supine



Left Lateral Prone



Lithotomy



Proctology

Brewer's unique approach to product development is centered on researching the needs of healthcare professionals and their patients. Before a design begins, we spend countless hours talking to doctors, nurses, and patients to understand what is most important to them. The result is a superior product that delivers the exact features needed to transform the procedure room into a safe, efficient, and comfortable environment. We invite you to discover the difference...

BrewerAssist

POWER PROCEDURE TABLE Infinitely adjustable for more possibilities.

Durability

- 3-year warranty
- Electric power system

Safety

- 19" low height
- Safety extra-low voltage
- Safety release actuators
- 450 lb. weight capacity*

Function

- 19" lowest height
- 40" maximum height
- 11.6" extendable footrest
- Adjustable three-point pivot headrest
- Intuitive foot and hand controls
- Easy to clean seamless upholstery
- 18" and 21" paper dispenser

Patient Comfort

- Easy lateral wheelchair transfer
- Infinitely adjustable headrest
- Extendable leg rest
- Multi-density, bolstered support upholstery
- Adjustable ergonomic stirrups

Aesthetics

- 18 standard upholstery color options
- Custom colors available

Weight Capacity

- 450 lbs.*

Table Dimensions

- Height: 19" - 40"
- Length: 69" - with footrest and headrest extends to 90"
- Width: 28"
- Leg rest: 16.5" x 11.25"
- Headrest: 13" x 11"
- Tilt range / Trendelenburg: to 30°

Paper Roll Holder

- Accepts 18" and 21" table paper rolls

Certifications

- UL listed, CUL listed

Electrical

- 115 VAC nominal, 60 Hz
- Safety extra-low voltage

Foot Control and Hand Pendant

- Foot control (standard)
- Hand pendant (optional)

#7000

accessories

- #99500 hand pendant
- #99501 side seat rails (Qty. 2)
- #99502 accessory clamp
- #99503 standard arm board
- #99504 deluxe articulating knee crutches (R&L)**
- #99505 stainless steel treatment pan
- #100253 urology drain pan system

*under normal use

**side seat rails and accessory clamp needed

All Brewer Assist Power Procedure Tables come with a 3-year warranty.

The Brewer Company reserves the right to improve or modify the design of its products consistent with new methods and techniques.

Call your local Brewer sales representative to learn about our complete product line.

The Brewer Company
13901 Main Street, Menomonee Falls, WI 53051
1.888.Brewer.1 www.brewercompany.com

Brewer
more possibilities