



**Brewer**  
more possibilities

**BASIC EXAM TABLE**

## The Brewer Basic Exam Table

### is clearly the:

- most reliable
  - most patient-friendly
  - most user-friendly
- entry exam table in the industry.

### See why Brewer means... more possibilities.

#### Durability

- 3-year warranty

#### Safety

- Sturdy, retractable patient step with mat
- 325 lb. weight capacity\*

#### Function

- Simple pneumatic back for easy patient positioning
- Easy to clean seamless upholstery
- Largest storage capacity available
- Front and side drawers
- Leg extension
- Plastic debris tray
- Pelvic tilt capability (optional)

#### Patient Comfort

- Sturdy step surface with mat
- Multi-density, bolstered support upholstery
- Adjustable ergonomic stirrups

#### Weight Capacity

- 325 lbs.\*

#### Table Dimensions

- Height: 32"
- Length: 57.5" extends to 71"
- Width: 26.5"
- Leg extension length: 15"

#### Aesthetics

- 18 standard upholstery color options
- Custom colors available



#### Paper Roll Holder

- Accepts 21" x 3.5" and 18" x 3.5"

#### Storage

- Front drawers: 19.7" x 19.2" x 4"
- Side drawers: 22.7" x 20.2" x 7"
- Side drawers available on right or left side

#### Certifications

- UL listed, CUL listed

#### Electrical Models

- 120 VAC, 60 Hz, 2 AMPS

#### Options

- Hospital-grade duplex electrical outlet available

#4000

#4001 with drawer heater and pelvic tilt

#### accessories

#99505 stainless steel treatment pan  
(12.75" x 10.38" x 2.5")

#100407 front drawer dividers

\*under normal use

### All Brewer Basic Exam Tables come with a 3-year warranty.

The Brewer Company reserves the right to improve or modify the design of its products consistent with new methods and techniques.

**Call your local Brewer sales representative to learn about our complete product line.**

The Brewer Company  
13901 Main Street, Menomonee Falls, WI 53051  
1.888.Brewer.1 www.brewercompany.com

**Brewer**  
more possibilities