

Alleviate™

Support System

The Science Behind Alleviate™.

There are as many as 300,000 musculoskeletal injuries among dental professionals each year associated with prolonged, awkward posture while performing procedures. While both irritating and painful, these injuries also account for lost time at work, a curtailment of normal activities, and for some, premature retirement.

With 60 years of manufacturing experience behind them, Brewer set out to provide better solutions to the painful problem. After much research and discussions with dentists and hygienists across the country, a support system was designed to help *alleviate* some of the musculoskeletal injuries within the dental industry.

When one or both arms are elevated away from the body, more weight is placed on the shoulders – adding pressure on the spine and muscles. Ideal posture includes shoulders relaxed with elbows close to the body. That's exactly what the Alleviate Support System is designed to do – relieve and reduce the occurrence of musculoskeletal pain and localized fatigue caused by excess weight on shoulders.

To order the Alleviate Support System or to schedule a demonstration, contact your local equipment dealer.



5" gel elbow rests provide comfort and support while absorbing vibration from instruments



patented ball and socket allows for free and natural movement



Alleviate™

Support System

tension control gives optimal flexibility and balance



versatile height adjustment of up to 7" enables a customized fit



9200 Series*



9100 Series



9000 Series



3300 Series

	height	seamless upholstery
9200BE	18.5"-24.5"	yes
9200BEST	18.5"-24.5"	-
9200BE-L	15.5"-20"	yes
9200BEST-L	15.5"-20"	-
9200BE-H	20.5"-26.5"	yes
9200BEST-H	20.5"-26.5"	-
9100BE	18.5"-24.5"	yes
9100BVE	18.5"-24.5"	yes
9100BE-L	15.5"-20"	yes
9100BVE-L	15.5"-20"	yes
9100BE-H	20.5"-26.5"	yes
9100BVE-H	20.5"-26.5"	yes

	height	seamless upholstery
9000BE	18.5"-24.5"	yes
9000BEST	18.5"-24.5"	-
9000BE-L	15.5"-20"	yes
9000BEST-L	15.5"-20"	-
9000BE-H	20.5"-26.5"	yes
9000BEST-H	20.5"-26.5"	-
3335BE	18"-25"	yes
3335BVE	18"-25"	yes

99444 - Alleviate (sold separately)

The Brewer Design Alleviate™ comes with a 5-year warranty.

Brewer Design reserves the right to improve or modify the design of its products consistent with new methods and techniques.

To order the Alleviate Support System or to schedule a demonstration, contact your local equipment dealer to learn about or complete product line.

Brewer Design
13901 Main Street, Menomonee Falls, WI 53051
1.800.558.8934 www.brewer-design.com

Brewer
DESIGN