

Clinically Researched Ergonomic Seating

Since 1947



Brewer Design, the most trusted name in dental seating, offers a variety of operator and hygiene seating options. A division of The Brewer Company, known for innovative and technologically advanced health care products, our product offers ergonomic support for your entire dental team.

Exclusively for Government Facilities, the Brewer program includes:

- 21-day lead time on 1-49 units
- 30-day lead time on 50+ units
- Dock-to-dock freight within the 48 contiguous states
- 5-Year complete repair/replace warranty

Quantity Discounts

- 10-20 units 5% off FSS Net Price
- 21-39 units 10% off FSS Net Price
- 40+ units 15% off FSS Net Price

Purchase Brewer products through Authorized distributors that fulfill your “set-aside” needs, including:

- Woman owned
- Service disabled veteran owned small business
- Accessibility across the U.S. for on-site demonstrations of Brewer stools.

**Access the entire Brewer FSS contract online at www.gsadvantage.gov
(contract #V797P - 3206M)**

Brewer Design Operatory Seating Selection and Specification Guide

9000 Series

Plush design with superior ergonomic features

- Cut out, saddle style seat measures 18" x 18"
- Large backrest, with lumbar support measures 18" x 18"
- Floating, synchronized seat and back tilt
- Standard height adjustment of 18.5" - 24.5" (Low and Higher heights available)
- Seamless upholstery in vacuum form materials only

Assistant Model

- Same features above
- Pneumatic height adjustment of 24" - 32"
- Ratcheted body support with adjustable foot ring



9100 Series

Our number one seller! Superior ergonomic support, designed with the smaller structured professional in mind. Great for smaller operatories.

- Smaller, saddle style seat measures 18" wide by 16" depth
- Cut out backrest, offers comfortable lumbar support
- Floating, synchronized seat and back tilt
- Standard height adjustment of 18.5" - 24.5" (Low and Higher heights available)
- Seamless upholstery with any material selection

Assistant Model

- Same features above
- Pneumatic height adjustment of 24" - 32"
- Ratcheted body support with adjustable foot ring



9200 Series

Complete ergonomic support and adjustability with standard size seat and smaller backrest.

- Saddle style seat measures 18" x 18"
- Cut out backrest, offers comfortable lumbar support
- Floating, synchronized seat and back tilt
- Standard height adjustment of 18.5" - 24.5" (Low and Higher heights available)
- Seamless upholstery in vacuum form materials only

Assistant Model

- Same features above
- Pneumatic height adjustment of 24" - 32"
- Ratcheted body support with adjustable foot ring



3300 Series

Traditional style seating offers ergonomic support while working long hours.

- Contoured, waterfall seat measures 18" x 18"
- Fully adjustable, contoured backrest measures 18" x 18"
- Synchronized seat and back tilt
- Standard height adjustment of 18"- 25" (Higher height available)
- Seamless upholstery with any material selection

Assistant Model

- Same upholstery as above
- Pneumatic height adjustment of 22"- 31"
- Ratcheted body support with adjustable foot ring
- Seat tilt model available



3100 Series

Traditional style seating offers basic ergonomic support while working long hours.

- Contoured, waterfall seat measures 18" x 18"
- Fully adjustable, contoured backrest measures 18" x 18"
- Dual adjustable, tilting backrest
- Standard height adjustment of 17.5"- 24" (Higher height available)
- Seamless upholstery with any material selection

Assistant Model

- Same upholstery as above without backrest
- Pneumatic height adjustment of 22"- 31"



2000 Series

Standard operatory seating offers traditional comfort.

- Comfortable round seat is 16" in diameter with 4" of high density foam
- Adjustable, contoured backrest measures 16" wide by 6" high
- Standard height adjustment of 18"- 24" (Higher height available)
- Seamless upholstery in vacuum form materials only

Assistant Model

- Same upholstery as above
- Available with or without backrest
- Pneumatic height adjustment of 20"- 26"
- 360° "Swing-a-Matic" body support on backless model (ratched on back model)
- Adjustable footrest standard



Alleviate Support System

Mobile, elbow support system minimizes stress on neck and shoulder muscles while working.

- 4" gel pad offers comfortable support
- Arms adjust 7" in height for custom fit
- Mobile support allows use while performing most procedures, unlike stationary arms
- Available on all 9000, 9100, 9200 and 3300 Series



Government Seating Catalog

Access the entire Brewer FSS contract online at
www.gsadvantage.gov
(contract #V797P - 3206M)

Brewer Design reserves the right to improve or modify the design of its products consistent with new methods and techniques.

Call your authorized Brewer Design dealer representative to learn about our complete line.

Brewer Design
13901 Main Street, Menomonee Falls, WI 53051
1.800.558.8934 www.brewer-design.com
Fax: 262.251.4391

Brewer
DESIGN
Division of The Brewer Company, LLC