




Brewer
more possibilities

IV POLES

IV POLES

SERIES Stable. Versatile. Reliable.




Aluminum Base Five-Leg

- Five-leg, 25" cast aluminum base
- Chrome plated steel tubing
- Dual wheel nylon casters
- Exclusive "No Loss" knob
- Height: 52.25" - 92.75"
- Load capacity: 45 lbs.
- Some assembly required

#11350
2 ram's horn hook

#11360
4 ram's horn hook



Heavy Base Four-Leg

- Tough four-leg, epoxy coated 22" steel base
- Chrome plated steel tubing
- Rubber wheel, ball bearing casters
- Exclusive "No Loss" knob
- Rubber bumpers on legs to protect walls
- Height: 51.75" - 93.5"
- Load capacity: 45 lbs.
- Some assembly required

#43403
2 ram's horn hook


#43406
4 ram's horn hook



Two-Hook Four-Leg

- Four-leg, epoxy coated 22" steel base
- Chrome plated steel tubing
- Dual wheel nylon casters
- Easy twist lock height adjustment
- Height: 51.25" - 91.75"
- Load capacity: 45 lbs.
- Some assembly required

#11300
2 ram's horn hook




Plastic Base Five-Leg

- Five-leg, 21" plastic base
- Chrome plated steel tubing
- Dual wheel nylon casters
- Exclusive "No Loss" knob
- Height: 51.75" - 94"
- Load capacity: 45 lbs.
- Some assembly required

#11320
2 ram's horn hook

#11321
4 ram's horn hook



Welded Hook Four-Leg

- Four-leg, epoxy coated 23" steel base
- Chrome plated steel tubing
- Rubber wheel, ball bearing casters
- Easy twist lock height adjustment
- Height: 47" - 81"
- Load capacity: 45 lbs.
- Some assembly required

#11290
2 ram's horn hook

All Brewer IV poles come with a 5-year warranty.

The Brewer Company reserves the right to improve or modify the design of its products consistent with new methods and techniques.

Call your local Brewer sales representative to learn about our complete product line.

The Brewer Company
13901 Main Street, Menomonee Falls, WI 53051
1.888.Brewer.1 www.brewercompany.com

Brewer
more possibilities