



Brewer
more possibilities

OVERBED TABLES

OVERBED TABLES

SERIES

Versatile. Convenient. Classic.



"H" Base Tables

- Easy spring loaded height adjustment
- Four point H-base with two locking dual wheel nylon casters
- Height: 28.25" - 43.25"
- Base: 15.25" x 26.38"; Top: 15" x 31"
- Weight capacity: 50 lbs.
- Some assembly required

#11600 "H" base
gray laminate top, gray epoxy finish on base

#11610 "H" base
walnut grain laminate top, chrome base



"H" Base Tilt-Top Tables

- Convenient tilt-top feature allows for easy reading and writing
- Lipped edge contains spills and keeps personal items from rolling off
- Easy spring loaded height adjustment
- Four point H-base with two locking dual wheel nylon casters
- Height: 29.5" - 46"
- Base: 15.25" x 26.38"; Top: 15" x 30"
- 0.25" lip
- Weight capacity: 25 lbs.
- Some assembly required

#11650 "H" base
gray laminate top, gray epoxy finish on base

#11660 "H" base
walnut grain laminate top, chrome base



"U" Base Tables

- "U" base to fit easily in restricted space areas
- Easy spring loaded height adjustment
- Four point U-base with two locking dual wheel nylon casters
- Height: 28.25" - 43.25"
- Base: 15.25" x 26.38"; Top: 15" x 31"
- Weight capacity: 50 lbs.
- Some assembly required

#11620 "U" base
walnut grain laminate top, chrome base

#11630 "U" base
gray laminate top, gray epoxy finish on base



"U" Base Tilt-Top Tables

- Convenient tilt-top feature allows for easy reading and writing
- Lipped edge contains spills and keeps personal items from rolling off
- "U" base to fit easily in restricted space areas
- Easy spring loaded height adjustment
- Four point U-base with two locking dual wheel nylon casters
- Height: 29.5" - 46"
- Base: 15.25" x 26.38"; Top: 15" x 30"
- 0.25" lip
- Weight capacity: 25 lbs.
- Some assembly required

#11670 "U" base
gray laminate top, gray epoxy finish on base

#11680 "U" base
walnut grain laminate top, chrome base



Low "H" Base Table

- Convenient table solution for low care beds
- Easy spring loaded height adjustment
- Four point H-base with two locking dual wheel nylon casters
- 1.75" casters and plastic glides included
- Height: Casters 19.75" - 26.75"; Glides 17.5" - 25.5"
- Base: 15.25" x 26.38"; Top: 15" x 31"
- Ground clearance: 1.8" with glides; 3.1" with casters
- Weight capacity: 50 lbs.
- Some assembly required

#11640 low "H" base
walnut grain laminate top, chrome base

All Brewer overbed tables come with a 5-year warranty.

The Brewer Company reserves the right to improve or modify the design of its products consistent with new methods and techniques.

Call your local Brewer sales representative to learn about our complete product line.

Brewer
more possibilities