

BREWER. CHOICE.

Brewer continues to expand our offering allowing you to streamline your ordering process with better products, delivered faster, at the best price. Choice is good.



Brewer Mobile Cart #63100

An economical option that allows 50% more storage with all the possibilities you'd expect from Brewer.

- Powder coated steel cabinet
- 3 small drawers, 1 large drawer
- Four black, locking casters
- Rubber bumpers prevent wall damage
- Contoured top covered in a bacteriological, chemical and fungal material
- Two handles integrated directly into the top surface
- Drawer dividers available (optional — item #100407)
- 5 year warranty
- 7 day lead time

1.888.Brewer.1
www.brewercompany.com

Brewer
more possibilities