

BrewerAccess

EXAM TABLE

Safety and efficiency for more possibilities.

#5000 Brewer Access Exam Table
#5001 with drawer heater and pelvic tilt

accessories

#99505 stainless steel treatment pan
(12.75" x 10.38" x 2.5")
#100407 front drawer dividers

Brewer
more possibilities

Durability

- 5-year warranty

Safety

- Large patient step
- 500 lb. weight capacity

Function

- Stable 4-leg patient step
- Patented pass-through side drawers
- Largest storage capacity available
- Additional paper storage under backrest
- Simple pneumatic backrest for easy patient positioning
- Easy to clean seamless upholstery
- Leg extension
- Plastic debris tray
- Pelvic tilt capability (optional)

Patient Comfort

- Largest step surface
- Multi-density, bolstered support upholstery
- Adjustable ergonomic stirrups

Aesthetics

- 18 standard upholstery color options
- 6 premium plush upholstery color options
- Custom colors available

Patented Versatility

- Pass-through side drawers can be opened either side
- Drawer stops [provided] allow drawer to open from only one side if desired
- Outlets can be reversed from right to left side

Weight Capacity

- 500 lbs. (under normal use)

Table Dimensions

- Height: 32"
- Length: 57.3" extends to 71"
- Width: 28"
- Leg extension length: 15"

Paper Roll Holder

- Accepts 21" x 3.5" and 18" x 3.5"
- Storage for up to 6 additional rolls under top

Storage

- Front drawers: 19.7" x 19.2" x 4"
- Side drawers: 22.7" x 20.2" x 7"
- Pass-through side drawers accessible from right or left side

Patient Step Dimensions

- Step width: 20.4"
- Step length: 17.7"
- Step height: 8"

Certifications

- UL listed, CSA listed

Electrical

- 120 VAC, 60 Hz, 10 AMPS
- Hospital-grade duplex electrical outlets