

BrewerAccess High-Low EXAM TABLE

Designed to take patients to a new height of care.

- #6000 pneumatic back high-low exam table
- #6001 pneumatic back with pelvic tilt, drawer heater
- #6500 power back high-low exam table
- #6501 power back with pelvic tilt, drawer heater

accessories

- #100407 front drawer dividers
- #98464 stainless steel treatment pan
- #100220 urology drain pan system
- #99504 deluxe articulating knee crutches
- #100230 deluxe articulating knee crutch clamp

Brewer
more possibilities

Durability

- 3-year warranty

Safety

- 18" low height
- Easy on / off for patients
- Safety grab bars for patient stability and positioning
- 450 lb. weight capacity
- Safety, extra low voltage
- Standard on / off switch

Function

- Patented pass-through side drawers
- Largest storage capacity available
- Pass-through work surface on side of table
- Backrest adjusts for easy patient positioning (pneumatic or power back available)
- Additional paper storage
- Return to chair feature (on power back model)
- Easy to clean, seamless upholstery
- Standard foot control

Patient Comfort

- Easy lateral wheelchair transfers
- Safety grab bars
- Multi-density support upholstery
- Adjustable ergonomic stirrups
- Toe board / leg extension / plastic debris tray
- Premium, plush upholstery available

Aesthetics

- 18 standard upholstery color options
- 6 premium, plush upholstery color options
- Custom colors available

Weight Capacity

- 450 lbs. (under normal use)

Table Dimensions

- Height: 18" - 37"
- Length: 57" extends to 72"
- Upholstery width: 28"
- Leg extension length: 15"

Covered Paper Roll Storage / Paper Roll Holder

- Storage capacity: 23" x 23" x 3.25"
- Storage for up to 6 additional paper rolls
- Accepts 21" x 3" and 18" x 3" paper rolls

Storage

- Front drawer: 19.7" x 19.2" x 4"
- Side drawers: 22.7" x 20.2" x 7"
- Pass-through side drawers accessible from right or left side of table
- Covered paper roll storage

Certifications

- UL listed, CUL listed

Options

- Power back / return to chair
- Premium plush upholstery
- 5° pelvic tilt / 60-Watt drawer heater in front top drawer