

BrewerAssist

POWER PROCEDURE TABLE

Infinitely adjustable for more possibilities.

#7000 Brewer Assist Power Procedure Table

accessories

#99500 hand pendant
#99501 side seat rails (Qty. 2)
#99502 accessory clamp
#99503 standard arm board
#99504 deluxe articulating knee crutches (R&L)*
#99505 stainless steel treatment pan
#100253 urology drain pan system

*side seat rails and accessory clamp needed

Durability

- 3-year warranty
- Electric power system

Safety

- 19" low height
- Safety extra-low voltage
- Safety release actuators
- 450 lb. weight capacity

Function

- 19" lowest height
- 40" maximum height
- 11.6" extendable footrest
- Adjustable three-point pivot headrest
- Intuitive foot and hand controls
- Easy to clean seamless upholstery
- 18" and 21" paper dispenser

Patient Comfort

- Easy lateral wheelchair transfer
- Infinitely adjustable headrest
- Extendable leg rest
- Multi-density, bolstered support upholstery
- Adjustable ergonomic stirrups

Aesthetics

- 18 standard upholstery color options
- 6 premium plush upholstery color options
- Custom colors available

Weight Capacity

- 450 lbs. (under normal use)

Table Dimensions

- Height: 19" - 40"
- Length: 69" - with footrest and headrest extends to 90"
- Width: 28"
- Leg rest: 16.5" x 11.25"
- Headrest: 13" x 11"
- Tilt range / Trendelenburg: to 30°

Paper Roll Holder

- Accepts 18" and 21" table paper rolls

Certifications

- UL listed, CUL listed

Electrical

- 115 VAC nominal, 60 Hz
- Safety extra-low voltage

Foot Control and Hand Pendant

- Foot control (standard)
- Hand pendant (optional)

Brewer
more possibilities